

[illegible]

cfishman@illumina.com

illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CSpPro, GenomeStudio, Genetic Energy, HiSeq, and HiScan are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

[illegible]

ILMN




Illumina Headquarters, San Diego, CA

illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CSpPro, GenomeStudio, Genetic Energy, HiSeq, and HiScan are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.



Historical Timeline (in brief)



Typical Genome Center in the Year 2006

- ▶ Large Sequencing Factories
- ▶ 100 to 150 Capillary Sequencers
- ▶ 5 to 10 Colony Picking Robots
- ▶ Dozens of PCR Machines
- ▶ Several Liquid Handling Robots
- ▶ Thousands of 384 well plates
- ▶ Dozens of lab personnel
- ▶ Multi-million dollar budgets



Broad, WashU, Sanger,
Baylor, BGI, Venter, JGI, etc.




Output of
3 to 5 Gb per month

Key Developments Over the Last 12 Months

Substantial enhancements on the Genome Analyzer

- >4x increase in output
- 2x increase in read length
- >2x the tags per run
- Further streamlined workflow



An Illumina Sequencer for Everyone

Next Generation
Sequencing made
accessible.
40GB/run



GA_{IIe}

Two proven
technologies. One
powerful platform.
50GB/run



HiScanSQ

Most widely adopted
NGS platform.
95GB/run



GA_{IIx}


Redefining the
trajectory of
sequencing.
200GB/run



HiSeq2000


Driving down the cost of human genome sequencing

Reagents price per 30x human genome




Publications

Cumulative Genome Analyzer Publications




GTATCATTAAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT
AGTAACACACTTCTGTAAAGCTTAAGATTAATGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT
CGTGGACAGTAAAGACACTTCTGTAAAGCTTAAGATTAATGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT
CTTGAAGCTTAAGATTAATGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT
CAACGATTAATGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT
ATTAAAGCTTAAGATTAATGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT
GTATCATTAAGATTAATGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGAGCTAAATTAAGGTACCATTAAGAGCTACCGGTGAGGACGAAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAT

Publications on Genome Analyzer



Simplest Sequencing Workflow

SIMPLIFIED SAMPLE PREP




cBot CLUSTER GENERATION



HiSeq 2000 SEQUENCING



DATA PROCESSING & ANALYSIS



Parallel sample processing

Automated cluster generation

Automated sequencing

Simple, efficient data analysis

HiSeq 2000

Redefining the trajectory of sequencing

HIGHEST OUTPUT

Initially capable of up to 200 Gb per run

FASTEST DATA RATE

~25 Gb/day

7-8 days for 2 x 100 bp

HIGHEST NUMBER OF READS

One billion single-end reads*

Two billion paired-end reads*

*Based on one billion clusters passing filter



Dual Surface Imaging


Cutting-edge imaging technology

TDI line-scanning technology, four CCDs

Fastest scanning and imaging method


Eliminates Mode Scrambler

Automated real-time auto-focus



Images clusters grown on both surfaces of flow cell


Huge gain in number of reads and sequence output



The Power of Line Scanning

Maximizing data rate

	Point Imaging	Area Imaging	Line Imaging
Stage & filter movement delays	+	-	+
Data transfer delays	+	-	+
Practical data acquisition limit	-	+	+
Data quality/background rejection	+	-	+




Four CCD Image Capture *Maximizes Throughput*

Cameras capture data almost 100% of time

Dual flow cells - alternation of imaging and chemistry

Simultaneous vs. sequential imaging



GTATCATTAAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGTAACTAAATTAAGGTACCATTAAGAGCTACCGGTGTAAGCAGCAAAAGATGATAACAGTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAC
AGTAACACAGTCTGTGTTAAAGCTTAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGTAACTAAATTAAGGTACCATTAAGAGCTACCGGTGTAAGCAGCAAAAGATGATAACAGTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAC
CGTTGATCATTAAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGTAACTAAATTAAGGTACCATTAAGAGCTACCGGTGTAAGCAGCAAAAGATGATAACAGTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAC
CAATCATTAAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGTAACTAAATTAAGGTACCATTAAGAGCTACCGGTGTAAGCAGCAAAAGATGATAACAGTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAC
ATTAAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGTAACTAAATTAAGGTACCATTAAGAGCTACCGGTGTAAGCAGCAAAAGATGATAACAGTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAC
GTATCATTAAAGATTAATTGATCCACTGATTCACGGTACCGTTACCGAAGGTATCAATTGAGTAACTAAATTAAGGTACCATTAAGAGCTACCGGTGTAAGCAGCAAAAGATGATAACAGTAACACACTTCTGTAAAGCTTAAGCGAAGGTATCATTAAATTAAC

HiSeq 2000 Dual Flow Cell Design

Instrument scalability and experimental flexibility

TWO INDEPENDENT FLOW CELLS

Simultaneously run applications that require different read lengths

Run in single or dual flow cell mode

SIMPLE FLOW CELL LOADING

Flow cells held by vacuum

No oil needed

LED switch ensures correct connection



HiSeq 2000

Plug-and-play reagents

PRE-CONFIGURED SEQUENCING REAGENTS


Only two minutes hands-on time

Up to 200 cycles per flow cell

Color & bar-coded for tracking

Temperature-controlled compartment

Integrated paired-end fluidics



HiSeq Control Software

Touch screen user interface

SIMPLE OPERATION

Step-by-step run set up



Enter read length, single, paired-read, indexing

Real time run

Reagents usage

Remote monitoring

Multimedia help



World's most adopted and proven next-gen sequencer

Internal runs have generated >95Gb

Version 5 reagents - 2 x 150 bp reads

95G Data Examples

Internal runs

Run	Ave error rate	Yield
A	2.02%	109Gb
B	1.74%	93.7Gb
C	1.21%	94.5Gb

Customer data

Run	Ave error rate	Yield
D*	1.63%	94.8Gb
E**	1.4%	97.8 Gb

**Part of early Access programme; extrapolated from 4 lanes*

***Broad Institute, courtesy of Sheila Fisher*

HiSeq 2000

Comparison with the Genome Analyzer




	HiSeq 2000 (at launch)	GA _{IIx} (at 50G)	GA _{IIx} (at 95G)
Gb per run	150-200	50	95
Gb per day	20-25	5	7
Cluster density in KClusters/mm ² **	260-350	490	620
Read length	2 x100	2 x100	2 x150
Available surface area (mm ²)*	2880	510	510

*GA_{IIx} with single surface, single FC, HiSeq 2000 with dual surface, dual FC

**Clusters passing filter

In one sequencing run you can...




New Sequencing Sample Prep Automation*

FOR GENOMIC DNA & mRNA-SEQ

Increases scalability

Reduces pipetting errors and variability



PRE-CONFIGURED PROTOCOLS

Accelerates adoption and transition times

No extensive development and verification needed

Serviced and supported by Illumina

*Expected availability End 2010

cBot

New automated cluster generation system

EASE OF USE

Pre-configured, plug-and-play reagents

Simple touch screen operation

Barcode scanner for reagent and sample tracking

Real-time run monitoring and self-diagnostics

AUTOMATION

Intervention-free clonal amplification in 4 hours


QUALITY

10-20% increase in cluster density

Less variability vs. cluster station



Sample preparation



1. Fragment template



	Nebulization	Covaris	Hydro Shear
Technology	Mechanical shearing	Adaptive Focused Acoustics	Hydrodynamic Shearing
Advantages	Low Cost Ease of Use	Low Sample Loss Tight Size Distrib. @ >1Kb	Enables Mate Pair analysis
Limitations	Large sample Loss Broad Size Distribution	Expensive, \$40,000 (S2) Broad Size Distrib. @ <1Kb	PreClearing of samples required. Assembly calibration required; assemblies can clog Low Throughput
Minimum Length of Starting Material	2 Kb	600bp	3 Kb
Maximum Product Fragment Length	5 Kb	4 Kb	20 Kb
Automatable	NO	YES (E210)	NO
Applications	ReSeq, PE	ReSeq, PE, Targeted ReSeq	Mate Pair

2. Repair Ends

- ▶ Converts overhangs into blunt ends and phosphorylates 5' end
- ▶ Reagents:
 - dNTP, T4 DNA pol, Klenow
 - Kinase/ATP (T4 PNK)
- ▶ Simple enzymatic reaction
- ▶ 30 minutes incubation


3. Adenylate Ends

- ▶ Adds 'A' base to the 3' end of the blunt, phosphorylated DNA fragments
- ▶ Prevents formation of adapters dimers and concatemers
- ▶ Reagents:
 - 1 mM dATP, Klenow exo (3' to 5' exo minus)


4. Link adapters

Same general template architecture regardless of application (gDNA, RNA, miRNA, ChIP-seq, exon pull-down, etc.)


Single Read



Paired Read




Indexed Paired Read



5. Gel purify

Size selection

- ▶ remove excess adapter and adapter dimers
- ▶ tighten the range of fragment sizes
- ▶ cut out a ~2 mm wide gel slice



Excision range

- ▶ 200 bp insert size target (+/- 10%) for 2×75 bp or shorter
- ▶ in practice, this translates to a 2 mm gel slice at ~300 bp to account for the length of the adapter sequences flanking the inserts
- ▶ 300 bp or greater insert size target for 2×100 bp or longer (unless you intentionally want to sequence overlapping read pairs)

6. Template enrichment /Final Adapter ligation

► Recommended number of cycles

Starting DNA (µg)	Volume of ligation product to PCR (µl)	Number of Cycle
5	1	10
1	5	10
0.5	10	10/12
0.1	10	12
0.05	all*	12

7. QC library


Quantify your completed library using


- ▶ NanoDrop (ng/μl)
- ▶ Agilent Bioanalyzer (size confirmation)

Bioanalyzer traces:


lane 1: 300bp library

lane 2: 600bp library







Cluster generation




The flow cell


- ▶ Everything except sample preparation is completed on the flow cell-Self Contained
 - ▶ Cluster Generation
 - ▶ Sequencing by Synthesis Reaction
 - ▶ Up to 96 samples
 - ▶ Closed system
- 



Cluster generation – hybridization and amplification




Cluster generation – hybridization and amplification




Cluster generation for single read runs

- ▶ Linearization & Blocking
- ▶ Denaturation & Primer Hybridization




Summary of cluster generation for single reads



How does this differ for paired end sequencing?


- The answer is in the flow cell oligos

SINGLE READ



Periodate Linearization


PAIRED END




Linearization 1 Enzyme

Linearization 2 Enzyme


Paired end cluster generation



Paired end cluster generation




Paired end linearization




So far, it's all the same as Single Read **except** we use an ENZYME (Linearization Enzyme 1) to perform the linearization instead of a CHEMICAL (periodate)

Paired end sequencing – read 1 and read 2




Sequencing



GTATCATTAAAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA
AGTAACACACTTCTGTAAAGCTTAAGAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA
CGTGGACAGAGTAAGACACTTCTGTAAAGCTTAAGAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA
CTTAAAGCTTAAGAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA
CAAGTAAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA
ATTAAAGCTTAAGAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA
GTATCATTAAGATTAATTGATCCACTGATTCACGGTACCGTAAAGGAGGTATCAATTGAGAGCTAAATATTAAGGTACCATTAAGAGGTACCGGTGAGGACGGAAGAGATGATAACATTAACACACTTCTGTAAAGCTTAAGCGGAAGGTATCATTAAAGATTA

NexGen Sequencing vs Sanger

Sanger sequencing



one region or PCR product
(read length 350 - 850 bp)

Illumina sequencing





1 cycle = 1 bp
36 cycles = 36 bp

whole genome
(read length 18 - 100 bp)

reference genome


sample 1
sample 2
sample 3
sample 4
sample 5
sample 6
sample 7
sample 8
sample 9
sample 10




48




Genome Analyzer IIX - Imaging compartment







GA IIx - Imaging Compartment




GA IIx - Imaging Set Up




Genome Analyzer Iix

SCS - *the data collection software*




GAIIx - Paired end sequencing




Sequencing by synthesis


1. Incorporation



Sequencing by synthesis



Sequencing by synthesis



1. Incorporation
2. Scan
3. Cleavage

Genome Analyzer_{IIx} Performance Specifications

Performance Parameters

- 27.5 – 33 Gb of high quality data / run
- 2.7 – 3.3 Gb / day
- 276-338M reads per paired-end run
- 2 x 100 bp supported read length
- Raw Accuracy:
 - ≥ 98% (2 x 100)
 - ≥ 99% (2 x 50)
- Run Time:
 - 2 x 100 bp in 9.5 days
 - 2 x 50 bp in 5 days
 - 1 x 35 bp in 2 days



The ENCODE Project

- a project with a vision to identify functional genomic elements




Fig. 1. Functional genomic elements being identified by the ENCODE pilot phase. The indicated methods are being used to identify different types of functional elements in the human genome.

www.sciencemag.org SCIENCE VOL 306 22 OCTOBER 2004


CLiP-Seq



Analysis

© 2010 Illumina, Inc. All rights reserved.


Data Analysis Overview



Software

- ▶ SCS 2.6 (**S**equencing **C**ontrol **S**tudio) → *GA PC*
- ▶ RTA 1.6 (Integrated **R**eal **T**ime **A**nalysis) → *GA PC*
- ▶ OLB 1.6 (**O**ff-Line **B**asecaller) → *Server*
- ▶ CASAVA 1.6 (**C**onsensus **A**ssessment of **S**equence **A**nd **V**ariation) → *Server*
 - *GERALD* (**G**eneration of **R**ecursive **A**nalyses **L**inked by **D**ependency)
 - ELAND v2 (**E**fficient **L**arge-Scale **A**lignment of **N**ucleotide **D**atabases)
 - PhageAlign

Sequencing Control Software v2.6



Real Time Analysis-RTA

▶ About RTA

- State Machine
- Output of every process is a file that acts as the trigger for the next step
- Multi-threaded
- Works in the background

▶ Input:

- Image files (.tiff)

▶ Output:

- qSeq.txt files

▶ Provides Real-Time metrics through graphs

▶ Transfers data from GA computer to server



X

X

How many images?

Real Time Analysis (RTA)


► Firecrest:

- Image Analysis/Template generation or cluster positioning
- Registration
 - Aligning template cluster positions with image
- Intensity Extraction
- Color Matrix (Cross Talk adjustment)
- Phasing/Prephasing

► Bustard:

- Basecalling
- qSeq file generation

Real Time Analysis-RTA



- ▶ Qseq.txt
- ▶ Tab-delimited: easy to parse, easy to import into databases
- ▶ Split files per read on a read pair / multiple read run

OLB 1.6

- “Offline Basecaller” (OLB)
- OLB is simply the old “Pipeline” software, but without the Gerald module
- OLB contains image analysis (i.e., “Firecrest”) and basecalling (“Bustard”)
- OLB is Linux-based
- All OLB algorithms identical to RTA 1.6 algorithms (equivalence directly tested vs. RTA)
- OLB will be released as a “tarball” (compress Linux .tar file)
- OLB Users Guide a separate downloadable .pdf file via iCom
- DVD’s no longer provided (not for PL, OLB or SCS)


CASAVA 1.6 High Level Features Summary

- ▶ GERALD and ELANDv2 Aligner
 - Supports gapped and multiseed alignment.
 - Use of Custom quality calibration tables is supported
 - Workflow simplification with the introduction of static make files
- ▶ Demultiplexing
 - Support for indexed runs (RTA via Demultiplex.pl and multiplexedGERALD.pl)
 - Sample Sheet - simple .csv file implementation
- ▶ Updated Allele Caller
 - Significantly faster (approximately 20-fold)
 - Correctly handles reads with gapped alignments and 100bp reads
- ▶ Updated SNP Caller
 - SNP caller updated to increase heterozygous SNP sensitivity.

CASAVA 1.6 High Level Features Summary

- ▶ GROUPER – Indel Detector
 - Uses singleton/shadow read pairs to detect indels.
- ▶ RNA SEQ improvements
 - ReadBases (counting by number of bases) is now the default counting method
 - Normalization of RNA counts for exons and genes now based on RPKM (reads per kb per million aligned reads)
- ▶ Performance Improvements
 - Lower memory footprint and bandwidth
 - Removed chromosome name parsing
 - Chromosomes will be designated by the fasta reference filename

Overall Analysis Workflow



Genome Studio



DNA Sequencing Module (DS)

- Software tools can be analyzed to discover and confirm SNPs and chromosomal breakpoint regions. Visualization tools display consensus reads in the reassembled genome and indicate SNPs with colored letters

ChIP Sequencing Module (CS)


- Create global binding site maps of DNA-associated proteins. Differential binding levels between experimental groups can be identified by comparing sequences, regions, and peaks in table or chromosome views.

RNA Sequencing Module (RS)

- This module performs digital gene expression profiling by aggregating data from the Genome Analyzer Pipeline Software to count the abundance of reads falling within specific exons, genes, and splice junctions. The data are then graphically displayed as tables or plots with GenomeStudio Software

Genome Studio Viewer

- **An integrated genome browser**
- **Fast and easy data import**
- **Standard stacked alignment plots**
- **A common platform for sequencing and microarray applications**





Genome Studio Module

Seq	Index	Sample Name	Chr Id	Position	A Bases	C Bases	G Bases	T Bases	Call	Bases Used	Total Bases	Call Score
TAACC	0	3fc_DNA_p...	1	45162	0	0	0	5	T	5	5	12.70
AGGGT	1	3fc_DNA_p...	1	53534	5	0	0	0	A	5	5	11.11
AACCC	2	3fc_DNA_p...	1	744366	4	0	0	0	A	4	4	10.10
ACCCA	3	3fc_DNA_p...	1	806588	0	0	5	0	G	5	5	12.60
GGGTT	4	3fc_DNA_p...	1	914761	4	0	0	0	A	4	4	10.60
CCCAA	5	3fc_DNA_p...	1	1055720	0	4	0	0	C	4	4	10.40
GTAGG	6	3fc_DNA_p...	1	1200215	0	0	5	0	G	5	5	11.06
CTAAC	7	3fc_DNA_p...	1	1200223	0	0	0	5	T	5	5	10.86
CTAAC	8	3fc_DNA_p...	1	1572479	0	4	0	1	C	5	5	10.50
AGGGT	9	3fc_DNA_p...	1	1718547	0	0	5	0	G	5	5	12.90
AACCC	10	3fc_DNA_p...	1	1771205	5	0	0	0	A	5	5	10.65
GGTTA	11	3fc_DNA_p...	1	2010974	0	0	0	4	T	4	4	10.10
AGGGT	12	3fc_DNA_p...	1	2021217	0	0	0	4	T	4	4	10.60
CTCAA	13	3fc_DNA_p...	1	2128539	0	0	4	0	G	4	4	10.40
CCCTA	14	3fc_DNA_p...	1	2129343	0	0	0	5	T	5	5	11.41
GTTAG	15	3fc_DNA_p...	1	2273702	0	0	0	4	T	4	4	10.40
TAGGG	16	3fc_DNA_p...	1	2522484	0	0	0	4	T	4	4	10.20
CTCTA	17	3fc_DNA_p...	1	2699065	0	1	0	4	T	5	5	10.40
CTAAC	18	3fc_DNA_p...	1	2772900	0	5	0	0	C	5	5	11.81
GTTAG	19	3fc_DNA_p...	1	2777567	0	0	4	0	G	4	4	10.40
CTAAC	20	3fc_DNA_p...	1	2874398	4	0	0	0	A	4	4	10.60
GGTTA	21	3fc_DNA_p...	1	2967799	4	0	0	0	A	4	4	10.30
GTTAG	22	3fc_DNA_p...	1	3013819	0	0	4	0	G	4	4	10.50

Illumina evaluated and third Party Tools



third-party genome analyzer data analysis tools

The evolution of next-generation sequencing technology will continue to require collaboration between researchers, thought leaders, and industry. Illumina encourages the open exchange of information and ideas, and we are committed to rapidly delivering solutions that will maximize your success.

Below you will find tools developed by Illumina sequencing customers that you can use for your Genome Analyzer data. These software applications* are designed to be used for a sequencing applications such as resequencing, ChIP-Seq, and Digital Gene Expression. Our customers have made these tools freely available to you.

Genome Alignment Browsers

- **Gbrowse - Genomic Browsing**
Generic Model Organism Database Project
<http://www.gmod.org/wiki/index.php/Gbrowse>
- **UCSC Browser - Genome browsing and comprehensive annotation**
Generic Model Organism Database Project
<http://genome.ucsc.edu/goldenPath/help/customTrack.html>
- **Staden Tools (GAP4, TGAP) - Alignment and Visualization for Small Data S**
James Bonifield (initially developed by Rodger Staden), Wellcome Trust Sanger In
<http://sourceforge.net/projects/staden/>

Alignment and Polymorphism Detection

- **BFAST – Blat-like Fast Accurate Search Tool**
Nils Homer, Stanley F. Nelson and Barry Merriman, University of California, Los Angeles
<http://genome.ucla.edu/bfast>
- **MAQ – Mapping and Assembly with Quality**
Hang Li, Sanger Centre
<http://maq.sourceforge.net/maq-man.shtml>
- **Bowtie – An ultrafast memory-efficient short read aligner**
Ben Langmead and Cole Trapnell, Center for Bioinformatics and Computational Biology
University of Maryland
<http://bowtie-bio.sourceforge.net/>

Genomic Assembly

- **Velvet – De novo assembly of short reads**
Daniel Zerbino and Ewan Birney, EMBL-EBI
<http://www.ebi.ac.uk/~zerbino/velvet/>
- **SSAKE – Assembly of short reads**
Rene Warren, et al, British Columbia Cancer Agency
<http://bioinformatics.oxfordjournals.org/cgi/content/full/23/4/500>
- **Euler – Genomic Assembly**
Pavel Pevzner and Mark Chaisson, University of California, San Diego
<http://nbcrc.sdsc.edu/euler/>


ChIP Sequencing

- **ChIP-Seq Peak Finder**
Barbara Wold, Cal Tech and Rick Meyers, Stanford University
<http://woldlab.caltech.edu/html/software/>



Digital Gene Expression

- **Comparative Count Display**

How much can you do with just *one lane* of GA data?



3000X BRCA1+BRCA2, 12 samples per lane



Illumina Next-Gen Sequencing

Sequence on a scale never before possible

ENABLES DEEP SEQUENCING OF MANY LARGE, COMPLEX GENOMES AND TRANSCRIPTOMES

FUEL MAJOR STUDIES INTO CANCER AND OTHER COMPLEX GENETIC DISEASES

DRIVE TRANSLATIONAL MEDICINE

MAKE ROUTINE WHOLE GENOME SEQUENCING POSSIBLE

Human

Crop

Livestock

Consumer Genomics

